

Anti-Prmt8 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nr22

Item Number: CM0019124



Introduction

High-quality anti-PRMT8 rabbit polyclonal antibody for Western blot and ELISA. Detects human, mouse, and rat PRMT8. Suitable for studies on protein arginine methylation, chromatin regulation, and brain-specific functions. SWISS: Q9NR22.

[Learn More](#)

Parameter	Value
Product Name	PRMT8 Rabbit pAb
Remarks/Alias	HRMT1L3; HRMT1L4; PRMT8
Species	Human
GeneID (Human)	56341
GeneID	56341
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-380 of human PRMT8 (NP_062828.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NR22
Protein Weight	45 kDa
Shipping	Ice bag