


Anti-Prmt9 Rabbit Polyclonal Antibody For Wb, Elisa - Q6P2P2

Item Number: CM0018133



Introduction

Rabbit polyclonal antibody targeting human PRMT9, suitable for WB and ELISA. A reliable tool for investigating alternative splicing and protein arginine methylation in human samples.

[Learn More](#)

Parameter	Value
Product Name	PRMT9 Rabbit pAb
Remarks/Alias	PRMT10; PRMT9
Species	Human
Gene ID (Human)	90826
Gene ID	90826
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human PRMT9 (NP_612373.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6P2P2
Protein Weight	95kDa
Shipping	Ice bag