


Anti-Pro Caspase-1 + P10 + P12 Rabbit Mab For Wb, Ip, Elisa - P29466

Item Number: CM0012746



Introduction

Detects human pro-Caspase-1 and its p10/p12 cleaved forms in WB, IP, and ELISA. Rabbit monoclonal antibody with reactivity across human, mouse, and rat. Research-grade tool for inflammation, pyroptosis, and immunology studies.

[Learn More](#)

Parameter	Value
Product Name	pro Caspase-1 + p10 + p12 Rabbit mAb
Remarks/Alias	ICE; P45; IL1BC; Cleaved-Caspase-1
Species	Human
GeneID (Human)	834
GeneID	834
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29466
Protein Weight	45kDa
Shipping	Ice bag