


Anti-Procr (Epcr) Rabbit Polyclonal Antibody For Wb And Elisa - Q9Unn8

Item Number: CM0025427



Introduction

Rabbit polyclonal antibody targeting human PROCR (endothelial protein C receptor, EPCR), suitable for Western blot and ELISA. Cross-reacts with mouse. Ideal for coagulation research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	PROCR Rabbit pAb
Remarks/Alias	CCCA; EPCR; CCD41; PROCR
Species	Human
GeneID (Human)	10544
GeneID	10544
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-210 of human PROCR (NP_006395.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UNN8
Protein Weight	27kDa
Shipping	Ice bag