


Anti-Procr Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Unn8

Item Number: CM0017702



Introduction

Rabbit polyclonal antibody targeting human PROCR (endothelial protein C receptor). Validated for WB, IHC-P, and ELISA; cross-reacts with mouse and rat. Suitable for studies of blood coagulation and endothelial function.

[Learn More](#)

Parameter	Value
Product Name	PROCR Rabbit pAb
Remarks/Alias	CCCA; EPCR; CCD41; PROCR
Species	Human
GeneID (human)	10544
GeneID	10544
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-209 of human PROCR (NP_006395.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UNN8
Protein Weight	27kDa
Shipping	Ice bag