

Anti-ProdH Rabbit Monoclonal Antibody For Wb, Ip, Ihc-P, Elisa - O43272

Item Number: CM0005843

Introduction

CamelBio offers anti-PRODH rabbit monoclonal antibody (O43272) validated for WB, IP, IHC-P, and ELISA, targeting human proline dehydrogenase 1. Suitable for proline metabolism and mitochondrial function studies.

[Learn More](#)

Parameter	Value
Product Name	PRODH Rabbit mAb
Remarks/Alias	POX; PIG6; HSPOX2; PRODH1; PRODH2; TP53I6
Species	Human
GeneID (Human)	5625
GeneID	5625
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O43272
Protein Weight	68kDa
Shipping	Ice bag