

Anti-Proenkephalin-A Monoclonal Antibody For Wb, If-P, Elisa - P01210

Item Number: CM0006783



Introduction

Rabbit monoclonal antibody against human proenkephalin-A (PENK). Validated for WB, IF-P, and ELISA; cross-reacts with mouse and rat. Ideal for pain perception and opioid signaling research.

[Learn More](#)

Parameter	Value
Product Name	Proenkephalin-A Rabbit mAb
Remarks/Alias	PE; PENK-A; Proenkephalin-A
Species	Human
GeneID (Human)	5179
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-267 of human PENK (NP_001129162.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01210
Protein Weight	31kDa
Shipping	Ice bag