


Anti-Profilin1 Polyclonal Antibody For Wb And Elisa - P07737

Item Number: CM0014516



Introduction

Rabbit polyclonal antibody to human Profilin1 (SWISS: P07737), validated for WB and ELISA. Cross-reacts with mouse, rat, and monkey. Profilin1 regulates actin dynamics and cytoskeleton structure. Ideal for cytoskeletal research.

[Learn More](#)

Parameter	Specification
Product Name	Profilin1 Rabbit pAb
Remarks/Alias	ALS18; Profilin1
Species	Human
Gene ID (Human)	5216
Gene ID	5216
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human Profilin1 (NP_005013.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	P07737
Protein Weight	15kDa
Shipping	Ice bag