


Anti-Profilin1 Rabbit Monoclonal Antibody For Wb, Elisa - P07737

Item Number: CM0005744



Introduction

Profilin1 Rabbit Monoclonal Antibody for WB and ELISA, recognizing human, mouse, rat Profilin-1 (P07737). Validated for cytoskeleton and actin-binding studies, ideal for neurodegeneration and ALS research.

[Learn More](#)

Parameter	Value
Product Name	Profilin1 Rabbit mAb
Remarks/Alias	ALS18; Profilin1
Species	Human
Gene ID (Human)	5216
Gene ID	5216
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 61-140 of human PFN1 (P07737).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07737
Protein Weight	15kDa
Shipping	Ice bag