

Anti-Progastrin-Releasing Peptide Rabbit Polyclonal Antibody For Wb, Elisa - P07492

Item Number: CM0022497



Introduction

Rabbit polyclonal antibody targeting human progastrin-releasing peptide (GRP), predicted MW 16 kDa, for WB and ELISA; cross-reacts with mouse. Suitable for gastrointestinal hormone and itch signaling research.

[Learn More](#)

Parameter	Value
Product Name	Progastrin-releasing peptide Rabbit pAb
Remarks/Alias	BN; GRP-10; proGRP; preproGRP; Progastrin-releasing peptide
Species	Human
GeneID (Human)	2922
GeneID	2922
Immunogen	Recombinant protein of human progastrin-releasing peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P07492
Protein Weight	16kDa
Shipping	Ice bag