

Anti-Progesterone Receptor Rabbit Monoclonal Antibody For Wb, Elisa - P06401

Item Number: CM0003786



Introduction

Rabbit monoclonal antibody targeting human Progesterone Receptor (PR, NR3C3). Suitable for Western blotting (WB) and ELISA applications. Ideal for studying hormone signaling and reproductive tissue biology.

[Learn More](#)

Parameter	Value
Product Name	Progesterone Receptor Rabbit mAb
Remarks/Alias	PR; NR3C3; Progesterone Receptor
Species	Human
Gene ID (Human)	5241
Immunogen	Recombinant fusion protein corresponding to amino acids 600-700 of human Progesterone Receptor (P06401)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P06401
Protein Weight	99 kDa
Shipping	Ice bag