

Anti-Progesterone Receptor Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P06401

Item Number: CM0017355



Introduction

High-quality rabbit polyclonal antibody targeting human Progesterone Receptor (PR, NR3C3). Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for reproductive biology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Progesterone Receptor Rabbit pAb
Remarks/Alias	PR; NR3C3; Progesterone Receptor
Species	Human
Gene ID (Human)	5241
Gene ID	5241
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human Progesterone Receptor (NP_000917.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06401
Protein Weight	99kDa
Shipping	Ice bag