

Anti-Progranulin (Granulin) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P28799

Item Number: CM0018065



Introduction

Rabbit polyclonal antibody against human Progranulin (Granulin) validated for WB, IF/ICC, and ELISA. Cross-reacts with human and mouse. Suitable for inflammation, wound healing, and lysosomal function studies.

[Learn More](#)

Parameter	Value
Product Name	Granulin Rabbit pAb
Remarks/Alias	GEP, GP88, PEPI, PGRN, CLN11, PCDGF, Granulin
Species	Human
Gene ID	2896
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-310 of human Granulin (NP_002078.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P28799
Protein Weight	64kDa
Shipping	Ice bag