

Anti-Prohibitin 2 (Phb2) Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q99623

Item Number: CM0007858



Introduction

Rabbit monoclonal antibody against human Prohibitin 2 (PHB2), validated for WB, IHC-P, IF/ICC, IP, and ELISA. Detects endogenous PHB2 in human, mouse, and rat samples. Ideal for studying mitochondrial function, transcriptional regulation, and cell signaling.

[Learn More](#)

Parameter	Value
Product Name	Prohibitin 2 (PHB2) Rabbit mAb
Remarks/Alias	BAP, REA, p22, hBAP, Bap37, BCAP37, PNAS-141, Prohibitin 2 (PHB2)
Species	Human
Gene ID (Human)	11331
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human Prohibitin 2 (Prohibitin 2 (PHB2)) (Q99623).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99623
Protein Weight	33kDa
Shipping	Ice bag