

Anti-Prohibitin 2 (Phb2) Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q99623

Item Number: CM0027165



Introduction

Rabbit polyclonal antibody against human Prohibitin 2 (PHB2), validated for WB, IHC-P, IP, ELISA. Reacts with human, mouse, and rat samples. Ideal for mitochondrial function, transcriptional regulation, and cancer research studies.

[Learn More](#)

Parameter	Value
Product Name	Prohibitin 2 (PHB2) Rabbit pAb
Remarks/Alias	BAP; REA; p22; hBAP; Bap37; BCAP37; PNAS-141; Prohibitin 2 (PHB2)
Species	Human
Gene ID (Human)	11331
Gene ID	11331
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99623
Protein Weight	33kDa
Shipping	Ice bag