

Anti-Prohibitin Rabbit Monoclonal Antibody For Wb - P35232

Item Number: CM0002618



Introduction

Rabbit monoclonal antibody targeting human Prohibitin (PHB, P35232). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Detects a 30kDa protein. Suitable for mitochondrial, cell signaling, and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	Prohibitin Rabbit mAb
Remarks/Alias	PHB; HEL-215; HEL-S-54e; Prohibitin
Species	Human
GeneID	5245
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human Prohibitin (P35232)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35232
Protein Weight	30kDa
Shipping	Ice bag