

Anti-Prohibitin Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P35232

Item Number: CM0013003



Introduction

Rabbit polyclonal antibody against human Prohibitin (PHB), validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat, monkey. Suitable for mitochondrial and signaling research.

[Learn More](#)

Parameter	Value
Product Name	Prohibitin Rabbit pAb
Remarks/Alias	PHB; HEL-215; HEL-S-54e; Prohibitin
Species	Human
Gene ID (Human)	5245
Gene ID	5245
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 73-272 of human Prohibitin (NP_002625.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	P35232
Protein Weight	30kDa
Shipping	Ice bag