

Anti-Prohibitin Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P35232

Item Number: CM0014958



Introduction

Rabbit polyclonal antibody against human Prohibitin (PHB), suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studies on mitochondrial, nuclear, and plasma membrane functions of prohibitin.

[Learn More](#)

Parameter	Details
Product Name	Prohibitin Rabbit pAb
Remarks/Alias	PHB, PHB1, HEL-215, HEL-S-54e, Prohibitin
Species	Human
GeneID (human)	5245
GeneID	5245
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 73-272 of human Prohibitin (NP_002625.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35232
Protein Weight	30kDa
Shipping	Ice bag