

Anti-Prokineticin-2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Hc23

Item Number: CM0028357

Introduction

Rabbit polyclonal antibody against human Prokineticin-2, validated for WB, IHC-P, IF/ICC, and ELISA, with cross-reactivity to mouse and rat. Ideal for circadian rhythm and gastrointestinal studies.

[Learn More](#)

Parameter	Value
Product Name	PROK2 Rabbit pAb
Remarks/Alias	BV8, HH4, PK2, KAL4, MIT1, PROK2
Species	Human
GeneID (Human)	60675
GeneID	60675
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-129 of human PROK2 (NP_001119600.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HC23
Protein Weight	14kDa
Shipping	Ice bag