


Anti-Prolyl Hydroxylase Phd1 (Egln2) Monoclonal Antibody For Wb - Q96Ks0

Item Number: CM0009688



Introduction

Rabbit monoclonal antibody against human EGLN2 (PHD1/HIF-PH1), validated for Western blotting and ELISA. Cross-reacts with human, mouse, and rat. Suitable for studying cellular oxygen sensing, HIF degradation, and hypoxia-related pathways.

[Learn More](#)

Parameter	Value
Product Name	Prolyl hydroxylase PHD1 (EGLN2) Rabbit mAb
Remarks/Alias	EIT6; PHD1; EIT-6; HPH-1; HPH-3; HIFPH1; HIF-PH1; Prolyl hydroxylase PHD1 (EGLN2)
Species	Human
Gene ID (Human)	112398
Gene ID	112398
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Prolyl hydroxylase PHD1 (EGLN2) (Q96KS0)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96KS0
Protein Weight	44kDa
Shipping	Ice bag