


Anti-Prolyl Hydroxylase Phd1 (Egln2) Polyclonal Antibody For Wb, Ip, Elisa - Q96Ks0

Item Number: CM0024454



Introduction

Rabbit polyclonal antibody targeting human Prolyl hydroxylase EGLN2 (PHD1/HIF-PH1). Validated for WB, IP, ELISA. Cross-reacts with mouse and rat. Suitable for hypoxia and HIF pathway research.

[Learn More](#)

Parameter	Value
Product Name	Prolyl hydroxylase PHD1 (EGLN2) Rabbit pAb
Remarks/Alias	EIT6; PHD1; EIT-6; HPH-1; HPH-3; HIFPH1; HIF-PH1
Species	Human
GeneID	112398
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human Prolyl hydroxylase PHD1 (EGLN2) (NP_444274.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96KS0
Protein Weight	44kDa
Shipping	Ice bag