

Anti-Prop1 Rabbit Polyclonal Antibody For Wb And Elisa - 075360

Item Number: CM0021574

Introduction

Validate homeobox protein prophet of Pit-1 (PROP1) with this rabbit polyclonal antibody. Suitable for WB and ELISA, with cross-reactivity to mouse. Ideal for pituitary development and CPHD2 research.

[Learn More](#)

Parameter	Value
Product Name	PROP1 Rabbit pAb
Remarks/Alias	CPHD2; PROP-1; PROP1
Species	Human
GeneID (Human)	5626
GeneID	5626
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human PROP1 (NP_006252.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	075360
Protein Weight	25kDa
Shipping	Ice bag