

Anti-Proprotein Convertase 9 (Pcsk9) Rabbit Mab For Wb, Elisa - Q8Nbp7

Item Number: CM0005951



Introduction

Rabbit anti-PCSK9 monoclonal antibody for WB and ELISA. Recognizes human, mouse, and rat PCSK9. Suitable for studying cholesterol homeostasis, LDLR degradation, and cardiovascular disease. SWISS: Q8NBP7.

[Learn More](#)

Parameter	Value
Product Name	Proprotein Convertase 9(PCSK9) Rabbit mAb
Remarks/Alias	FH3; PC9; FHCL3; NARC1; LDLCQ1; NARC-1; HCHOLA3; Proprotein Convertase 9(PCSK9)
Species	Human
GeneID (Human)	255738
GeneID	255738
Immunogen	Synthetic peptide corresponding to amino acids 500-600 of human Proprotein Convertase 9(PCSK9) (Q8NBP7)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NBP7
Protein Weight	74 kDa
Shipping	Ice bag