

Anti-Proprotein Convertase 9 (Pcsk9) Rabbit Polyclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - Q8Nbp7

Item Number: CM0028127



Introduction

Rabbit polyclonal antibody against human Proprotein Convertase 9 (PCSK9). Suitable for WB, IP, IF/ICC, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Targets the 74 kDa cholesterol regulator PCSK9. Ideal for LDLR degradation and cholesterol metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Proprotein Convertase 9(PCSK9) Rabbit pAb
Remarks/Alias	FH3; PC9; FHCL3; NARC1; LDLCQ1; NARC-1; HCHOLA3; Proprotein Convertase 9(PCSK9)
Species	Human
Gene ID (Human)	255738
Gene ID	255738
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 590-692 of human Proprotein Convertase 9(PCSK9) (NP_777596.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NBP7
Protein Weight	74 kDa
Shipping	Ice bag