

Anti-Pros1 Rabbit Polyclonal Antibody For Wb, Elisa - P07225

Item Number: CM0014694



Introduction

High-quality rabbit polyclonal antibody against human PROS1 (Protein S). Validated for WB and ELISA with reactivity across human, mouse, and rat. Ideal for coagulation and thrombosis research.

[Learn More](#)

Parameter	Value
Product Name	PROS1 Rabbit pAb
Remarks/Alias	PSA; PROS; PS21; PS22; PS23; PS24; PS25; THPH5; THPH6; PROS1
Species	Human
GeneID (Human)	5627
GeneID	5627
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 405-676 of human PROS1 (NP_000304.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07225
Protein Weight	75kDa
Shipping	Ice bag