

Anti-Prostatic Acid Phosphatase (Acpp) Monoclonal Antibody For Wb And Elisa - P15309

Item Number: CM0009709



Introduction

Rabbit monoclonal antibody targeting human Prostatic Acid Phosphatase (ACPP), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for prostate cancer, phosphatase activity, and HIV SEVI fibril research.

[Learn More](#)

Parameter	Value
Product Name	Prostatic Acid Phosphatase (ACPP) Rabbit mAb
Remarks/Alias	ACPP; 5'-NT; ACP-3; TM-PAP; Prostatic Acid Phosphatase (ACPP)
Species	Human
GeneID (Human)	55
GeneID	55
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Prostatic Acid Phosphatase(PAP/ACPP) (P15309).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P15309
Protein Weight	45 kDa
Shipping	Ice bag