

Anti-Prostatic Acid Phosphatase (Acpp) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P15309

Item Number: CM0020104



Introduction

Rabbit polyclonal antibody against human Prostatic Acid Phosphatase (ACPP). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for cancer research and phosphatase studies.

[Learn More](#)

Parameter	Value
Product Name	Prostatic Acid Phosphatase (ACPP) Rabbit pAb
Remarks/Alias	ACPP; 5'-NT; ACP-3; TM-PAP; Prostatic Acid Phosphatase (ACPP)
Species	Human
Gene ID (Human)	55
Gene ID	55
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-386 of human PAP/Prostatic Acid Phosphatase (ACPP) (NP_001090.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15309
Protein Weight	45kDa
Shipping	Ice bag