

Anti-Prostatic Acid Phosphatase (Acpp) Rabbit Polymab® For Wb, Ihc-P, Elisa - P15309

Item Number: CM0003480



Introduction

Rabbit recombinant monoclonal antibody against human Prostatic Acid Phosphatase (ACPP), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for cancer and HIV studies.

[Learn More](#)

Parameter	Value
Product Name	Prostatic Acid Phosphatase (ACPP) Rabbit PolymAb®
Remarks/Alias	ACPP; 5'-NT; ACP-3; TM-PAP
Species	Human
GeneID (Human)	55
GeneID	55
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-386 of human Prostatic Acid Phosphatase (ACPP) (NP_001090.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15309
Protein Weight	45 kDa
Shipping	Ice bag