

Anti-Prostatic Acid Phosphatase (Psap) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P15309

Item Number: CM0003796



Introduction

Recombinant rabbit monoclonal antibody targeting human Prostatic Acid Phosphatase (PSAP/ACPP). Validated for WB, IHC-P, ELISA with cross-reactivity to human and rat. Suitable for prostate cancer biomarker research and HIV infection studies.

[Learn More](#)

Parameter	Value
Product Name	Prostatic Acid Phosphatase (PSAP) Rabbit mAb
Remarks/Alias	ACPP; 5'-NT; ACP-3; TM-PAP; Prostatic Acid Phosphatase (PSAP)
Species	Human
GeneID (Human)	55
GeneID	55
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 33-386 of human Prostatic Acid Phosphatase (ACPP) (NP_001090.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P15309
Protein Weight	45kDa
Shipping	Ice bag