


Anti-Proteasome 20S Lmp7/Psmb8 Monoclonal Antibody For Wb, Ihc-P, Elisa - P28062

Item Number: CM0005749



Introduction

Rabbit monoclonal antibody against proteasome subunit beta type-8 (LMP7/PSMB8) for WB, IHC-P, ELISA. Cross-reacts with Human, Mouse, Rat. Ideal for studying immunoproteasome function and antigen processing.

[Learn More](#)

Parameter	Value
Product Name	Proteasome 20S LMP7/PSMB8 Rabbit mAb
Remarks/Alias	JMP; ALDD; LMP7; NKJO; D6S216; PRAAS1; PSMB5i; RING10; D6S216E; Proteasome 20S LMP7/PSMB8
Species	Human
GeneID (Human)	5696
GeneID	5696
Immunogen	Recombinant protein of human Proteasome 20S LMP7/PSMB8
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28062
Protein Weight	30kDa
Shipping	Ice bag