


Anti-Proteinase K Rabbit Polyclonal Antibody - P06873

Item Number: CM0021463



Introduction

Rabbit polyclonal antibody against Proteinase K, validated for WB and ELISA. Specific to fungal Proteinase K (P06873), suitable for studying proteolytic activity and keratin hydrolysis.

[Learn More](#)

Parameter	Value
Product Name	Proteinase K Rabbit pAb
Species	Human
Immunogen	Recombinant protein of Proteinase K
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Fungi
SWISS	P06873
Shipping	Ice bag