

Anti-Prothrombin (F2) Rabbit Polyclonal Antibody For Wb And Elisa - P00734

Item Number: CM0017158



Introduction

High-quality rabbit polyclonal antibody targeting human prothrombin (F2). Validated for Western blot and ELISA, this antibody supports coagulation research, hemostasis studies, and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	F2 Rabbit pAb
Remarks/Alias	PT; THPH1; RPRGL2; F2
Species	Human
Gene ID (Human)	2147
Gene ID	2147
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 364-622 of human F2 (NP_000497.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P00734
Protein Weight	70kDa
Shipping	Ice bag