

Anti-Protoporphyrinogen Oxidase (Ppox) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P50336

Item Number: CM0026752



Introduction

Rabbit polyclonal antibody against human protoporphyrinogen oxidase (PPOX), validated for Western blot, immunohistochemistry (paraffin), and ELISA. Detects human and mouse PPOX. PPOX is a mitochondrial enzyme catalyzing a critical step in heme biosynthesis; mutations are linked to variegate porphyria.

[Learn More](#)

Parameter	Specification
Product Name	PPOX Rabbit pAb
Remarks/Alias	VP, PPO, V290M, PPOX
Species	Human
Gene ID (Human)	5498
Gene ID	5498
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human PPOX (NP_001116236.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P50336
Protein Weight	51kDa
Shipping	Ice bag