


Anti-Prox1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q92786

Item Number: CM0009640



Introduction

High-quality PROX1 Rabbit Monoclonal Antibody validated for WB, IHC-P, and ELISA applications. Reacts with human, mouse, and rat specimens. Target weight 83kDa. Ideal for developmental biology and circadian rhythm research.

[Learn More](#)

Parameter	Value
Product Name	PROX1 Rabbit mAb
Remarks/Alias	PROX1; prospero homeobox 1
Species	Human
GeneID (Human)	5629
GeneID	5629
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92786
Protein Weight	83kDa
Shipping	Ice bag