


Anti-Prox1 Rabbit Polyclonal Antibody For Wb, Elisa - Q92786

Item Number: CM0014302



Introduction

This PROX1 rabbit polyclonal antibody is validated for Western blot and ELISA, detecting human, mouse, and rat PROX1, and is ideal for developmental biology and circadian rhythm studies.

[Learn More](#)

Parameter	Value
Product Name	PROX1 Rabbit pAb
Remarks/Alias	PROX1
Species	Human
GeneID (Human)	5629
GeneID	5629
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 220-375 of human PROX1 (NP_002754.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92786
Protein Weight	83kDa
Shipping	Ice bag