

Anti-Prp19 Polyclonal Antibody For Wb, Elisa - Q9Ums4

Item Number: CM0028430



Introduction

Rabbit polyclonal antibody against human PRP19 (Pre-mRNA-processing factor 19), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for splicing and DNA repair studies.

[Learn More](#)

Parameter	Value
Product Name	PRP19 Rabbit pAb
Remarks/Alias	PSO4, SNEV, PRP19, UBOX4, hPSO4, NMP200
Species	Human
Gene ID (Human)	27339
Gene ID	27339
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human PRP19 (NP_055317.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UMS4
Protein Weight	55kDa
Shipping	Ice bag