

Anti-Prp19 Polyclonal Antibody For Wb - Q9Ums4

Item Number: CM0016538



Introduction

Rabbit polyclonal antibody targeting human PRP19 (PSO4/SNEV), validated for WB, IHC, IF/ICC, IP, and ELISA. Cross-reacts with human, mouse, rat. Ideal for studying pre-mRNA splicing and DNA damage repair mechanisms.

[Learn More](#)

Parameter	Value
Product Name	PRP19 Rabbit pAb
Remarks/Alias	PSO4; SNEV; PRP19; UBOX4; hPSO4; NMP200
Species	Human
GeneID (Human)	27339
GeneID	27339
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 127-416 of human PRP19 (NP_055317.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UMS4
Protein Weight	55kDa
Shipping	Ice bag