

Anti-Prp19 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Ums4

Item Number: CM0006247



Introduction

Recombinant rabbit monoclonal antibody targeting human PRPF19 (PRP19/SNEV). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for investigating pre-mRNA splicing, DNA repair, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	PRP19 Rabbit mAb
Remarks/Alias	PSO4; SNEV; PRP19; UBOX4; hPSO4; NMP200
Species	Human
GeneID (Human)	27339
GeneID	27339
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-504 of human PRPF19 (Q9UMS4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UMS4
Protein Weight	55kDa
Shipping	Ice bag