


Anti-Prp6 Rabbit Polyclonal Antibody For Wb, Elisa - O94906

Item Number: CM0023222



Introduction

Anti-PRP6 rabbit polyclonal antibody validated for WB and ELISA with human reactivity. Key tool for studying pre-mRNA splicing, androgen receptor signaling, and retinitis pigmentosa mechanisms. - O94906

[Learn More](#)

Parameter	Value
Product Name	PRPF6 Rabbit pAb
Remarks/Alias	TOM; ANT1; Prp6; RP60; ANT-1; hPrp6; U5-102K; C20orf14; SNRNP102; PRPF6
Species	Human
GeneID (Human)	24148
GeneID	24148
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O94906
Protein Weight	107kDa
Shipping	Ice bag