


Anti-Prpf3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - O43395

Item Number: CM0028849



Introduction

High-quality rabbit polyclonal antibody detecting PRPF3, a key spliceosome component. Validated in WB, IHC-P, IF/ICC, IP, ELISA with cross-reactivity to human, mouse, rat. Ideal for splicing and retinitis pigmentosa research.

[Learn More](#)

Parameter	Value
Product Name	PRPF3 Rabbit pAb
Remarks/Alias	PRP3, RP18, HPRP3, Prp3p, HPRP3P, SNRNP90, PRPF3
Species	Human
GeneID (Human)	9129
GeneID	9129
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human PRPF3 (NP_004689.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43395
Protein Weight	78kDa
Shipping	Ice bag