

Anti-Prpf4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O43172

Item Number: CM0028133



Introduction

Rabbit polyclonal antibody targeting PRPF4, a U4/U6 snRNP protein crucial for spliceosome assembly and pre-mRNA splicing. Validated for WB, IHC-P, and ELISA with human, mouse, rat reactivity. For research use.

[Learn More](#)

Parameter	Value
Product Name	PRPF4 Rabbit pAb
Remarks/Alias	PRP4, RP70, HPRP4, Prp4p, HPRP4P, SNRNP60, PRPF4
Species	Human
GeneID	9128
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human PRPF4 (NP_004688.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43172
Protein Weight	58kDa
Shipping	Ice bag