


Anti-Prpf4B Rabbit Polyclonal Antibody For Wb And Elisa - Q13523

Item Number: CM0014763



Introduction

PRPF4B Rabbit Polyclonal Antibody validated for WB and ELISA. Recognizes human, mouse, rat PRPF4B, a kinase involved in spliceosomal assembly, mitosis, and Hippo signaling. Ideal for chromatin and splicing research.

[Learn More](#)

Parameter	Value
Product Name	PRPF4B Rabbit pAb
Remarks/Alias	PR4H; PRP4; PRP4H; PRP4K; dJ1013A10.1; PRPF4B
Species	Human
Gene ID (Human)	8899
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 550-760 of human PRPF4B (NP_003904.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13523
Protein Weight	117kDa
Shipping	Ice bag