

Anti-Prpf8 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q6P2Q9

Item Number: CM0008188



Introduction

CamelBio provides PRPF8 Rabbit mAb (Q6P2Q9) validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat PRPF8, a key splicing factor. Recombinant immunogen ensures specificity. Ideal for pre-mRNA splicing research.

[Learn More](#)

Parameter	Value
Product Name	PRPF8 Rabbit mAb
Remarks/Alias	PRP8; RP13; HPRP8; PRPC8; SNRNP220; PRPF8
Species	Human
GeneID (Human)	10594
GeneID	10594
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2200-2335 of human PRPF8 (Q6P2Q9)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6P2Q9
Protein Weight	274kDa
Shipping	Ice bag