

Anti-Prph Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P41219

Item Number: CM0025794

Introduction

Rabbit polyclonal antibody against human PRPH (Peripherin), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for neuronal intermediate filament research and neurodegeneration studies.

[Learn More](#)

Parameter	Value
Product Name	PRPH Rabbit pAb
Remarks/Alias	NEF4; PRPH1; PRPH
Species	Human
Gene ID	5630 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human PRPH (NP_006253.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41219
Protein Weight	54kDa
Shipping	Ice bag