


Anti-Prps2 Polyclonal Antibody For Wb, Elisa - P11908

Item Number: CM0018556



Introduction

High-quality rabbit polyclonal antibody targeting human PRPS2 (ribose-phosphate pyrophosphokinase 2), suitable for WB and ELISA. Cross-reacts with rat. Key tool for nucleotide biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	PRPS2 Rabbit pAb
Remarks/Alias	PRSI; PRPS2
Species	Human
Gene ID (Human)	5634
Gene ID	5634
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human PRPS2 (NP_002756.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P11908
Protein Weight	35kDa
Shipping	Ice bag