


Anti-Prps2 Polyclonal Antibody For Wb, Elisa - P11908

Item Number: CM0013479



Introduction

Rabbit polyclonal antibody targeting human PRPS2 (ribose-phosphate pyrophosphokinase 2), validated for WB and ELISA. Detects human, mouse, and rat PRPS2. Ideal for nucleotide biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	PRPS2 Rabbit pAb
Remarks/Alias	PRSII; PRPS2
Species	Human
GeneID (Human)	5634
GeneID	5634
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human PRPS2 (NP_002756.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11908
Protein Weight	35kDa
Shipping	Ice bag