


Anti-Prpsap2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O60256

Item Number: CM0017301



Introduction

Rabbit polyclonal antibody against human PRPSAP2 (PAP41), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Detects a 41 kDa band. Suitable for nucleotide biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	PRPSAP2 Rabbit pAb
Remarks/Alias	PAP41; PRPSAP2
Species	Human
GeneID (Human)	5636
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 211-252 of human PRPSAP2 (NP_001340030.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60256
Protein Weight	41kDa
Shipping	Ice bag