

Anti-Prr16 Rabbit Polyclonal Antibody For Wb, Elisa - Q569H4

Item Number: CM0023122

Introduction

Rabbit polyclonal antibody targeting human PRR16/Largen, validated for WB and ELISA, with mouse cross-reactivity. Regulates cell size and mitochondrial function; ideal for cell biology and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	PRR16 Rabbit pAb
Remarks/Alias	DSC54; LARGEN; PRR16
Species	Human
GenelD (human)	51334
GenelD	51334
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q569H4
Protein Weight	33kDa
Shipping	Ice bag