

Anti-Prr5 Polyclonal Antibody For Wb, If/Icc, Elisa - P85299

Item Number: CM0013249

Introduction

Rabbit polyclonal antibody against PRR5 (Protor-1), an mTORC2 subunit and tumor suppressor, validated for WB, IF/ICC, and ELISA in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	PRR5 Rabbit pAb
Remarks/Alias	PP610; PROTOR1; PROTOR-1; FLJ20185k; PRR5
Species	Human
GeneID (human)	55615
GeneID	55615
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 190-388 of human PRR5 (NP_851850.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P85299
Protein Weight	43kDa
Shipping	Ice bag