

Anti-Prrx1 Rabbit Polyclonal Antibody For Wb And Elisa - P54821

Item Number: CM0015929

Introduction

Rabbit polyclonal antibody targeting human PRRX1 (PMX1, PRX1). Validated for Western blot (WB) and ELISA; cross-reacts with human and mouse. Suitable for research on transcriptional regulation, myofibroblast differentiation, and wound healing.

[Learn More](#)

Parameter	Value
Product Name	PRRX1 Rabbit pAb
Remarks/Alias	PMX1; PRX1; AGOTC; PHOX1; PRX-1; PRRX1
Species	Human
GeneID (Human)	5396
GeneID	5396
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human PRRX1 (NP_073207.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P54821
Protein Weight	27kDa
Shipping	Ice bag